

J4 mini head dental handpiece for children



Technical parameters


Rotate speed	Noise	Power	Water flow	Air pressure	Diameter of head (mini head)	Chucking type
350000-500000rpm	< 65dB	14W	90-110ml/min	0.23-0.25Mpa	8.9mm	Push button

Spare parts

Bearing NO: 301060	Spindle NO: 301199	Cartridge NO: 202151	Rear cover NO: 200181
			



Every dentist deserves to have a piece of mini head dental hanpiece for children — — J4



Guangdong JINME Medical Technology Co., Ltd.
 Address: A15, New Lighting Source Industrial Base, Foshan, China
 Tel: +86-757-88777979
<http://www.jinmedental.cn>
 E-mail: info@jinmedental.com



Distributor

J4 mini head dental handpiece for children

J4

J4 mini head dental handpiece for children

Mini head provides more flexible operation

8.9mm extraordinary miniature head leaves much more flexible and bigger operational space in children's oral cavity.



J4 mini head dental handpiece for children



Standard dental handpiece



German bearing is adopted

Cooperate with the advanced German supplier on miniature bearing so as to provides the gaurantee for the supper-long life span of the handpiece.



Equipped with professional short shank diamond burs for children's teeth

Make the dentist's treatment more convenient and precise with professional short shank diamond bur, it avoids the influence of children's limited intraoral operational space.